

Micro Sensor Background Measurement Frequency “Calibration”




neospectra
si-WARE
SYSTEMS

Background Measurement Need

- **Absorbance is measured by ratioing Sample PSD to Background PSD, so the instrument spectral response is compensated**
- **Frequent measurement of background is recommended to avoid effect of instrument response drift**
- **Scanner response drift**
 - **Self-Heating (Warm-up) Effect**
 - No environmental changes (Ambient Temperature, Humidity, ...), variations will come from self heating/warm-up of the scanner (Sensor + Other electronics in scanner)
 - It depends on light source usage frequency and electronics operation inside the scanner (always on or off)
 - **Environmental Changes (Ambient Temperature, Humidity)**
 - Self-correction (internal sensor calibration) should be applied and new background is needed
 - Highly unlikely in case of a salon with controlled temperature
- **Can be solved using some DSP in chemometrics (SNV or baseline correction algorithms)**

Self-Heating (Warm-up) Effect Overview

Example in Certain Conditions



- **Continuous measurements**
- **Ratioing measurements to first measurement**
- **Baseline correction (SNV simplest) can solve this issue**

Characterization of self-heating

- Characterization of baseline drifts considering frequency of usage (how many scans per hour ? Separation between measurements ?)
- Decide on: start-up initialization time needed

Effect of self-heating on prediction error

- Multiple measurements of a standard hair sample during self-heating period
- Multiple measurements of a standard hair sample across a day with different time separation between measurements
- Decide if we need for
 - Start-up time initialization time needed
 - Working on baseline correction algorithm

**Thank
you**



neospectra
— **si-ware**
SYSTEMS